



# SprintRay ProWash 8

Multistage Automated Washing System Designed for Digital Dentistry

## Technical Specifications April 2023

<b>Washing Technology</b>	Patented Mechanical Jetting Technology Automatic 2-stage wash and dry cycle Precision 10,000 RPM propeller	
<b>Solvent Management</b>	Visual representation of solvent usage and precalibrated Hydrometers Capacitive solvent level Sensor	
<b>Wash Volume/Envelop</b>	<b>Platform Envelope</b> 7 x 5 x 3 in. 17 x 12.7 x 7.6 cm	<b>Basket Envelope</b> 7 x 5 x 5 in. 17 x 12.7 x 12.7 cm
<b>Solvent Compatibility</b>	Isopropyl alcohol 91% Isopropyl alcohol 99%	
<b>Solvent Capacity</b>	Min 1L per reservoir Max 2L per reservoir	
<b>User Controls</b>	5 inch color touch-screen Predefined wash/dry preset User programmable wash/dry option	
<b>Connectivity</b>	Dual band 2.4/5 GHz Wi-Fi chipset Local network via ethernet cable USB-C port	
<b>Safety Sensors</b>	Automatic start/stop sensors Platform vs basket detection Fill level sensor	
<b>Unit Dimensions (W x D x H)</b>	12 x 12 x 18 in. 30 x 30 x 46 cm	
<b>Shipping Weight &amp; Dimensions</b>	14 x 14 x 20 in. 36 x 36 x 50 cm 45 lbs. / 20.4 kg	
<b>Power Requirement</b>	100-240 V. (via DC adapter)	

